

Molecular Pathology Laboratory Network, Inc.

Client Supply Request Form
Form Version March 2005

Client Account #

00000

Client Information

Request Date: _____

Requested By: _____

Phone Number: _____

Fax Number: _____

Please complete this form and fax to Molecular Pathology Laboratory Network, Inc. at 865-273-1128 when supplies are needed. To help us better serve you, please complete all contact information and check stock on all items before submitting this request. For questions or comments contact our Materials Management Specialist at 865-380-9746 (local) or toll-free at 1-800-932-2943, extension 167.

Supplies	Qty.	Comments
Small box (styrofoam with grey foam insert)- 6.5" x 5" x 8" *		
Medium flat cardboard shipping box (Item # 1162)- 11" x 6.5" x 2.75"		
Large cardboard box (dry ice)- 11" x 9" x 10" **		
DHL/ Airborne Express Preprinted Airbills		
DHL/ Airborne Express Lab Shipping Packs		
FedEx Preprinted Airbills		
FedEx Airbill Pouches		
FedEx Clinical Packs (formerly Diagnostic Packs)		
MPLN Test Order Requisitions		
Freezer Gel Packs		
Biohazard Specimen Transport Bags		
Pour-Off Tubes , 5ml screw cap (for diagnostic specimens)		
Digene HPV DNA Cervical Swab Collection Kit		
MPLN Pathogen Collection Kit (swab with transport medium)		
PPT Specimen Collection Tubes (HIV specimens only)		
RPMI Tissue Transport Medium- Standard 10ml Tubes		
RPMI Tissue Transport Medium- Small (FNA) 3ml Tubes		
Other (Please Specify)		

* Styrofoam shippers measuring approx. 6.5" x 5" x 8" with white paperboard outer sleeve. Please order in multiples of six (6) each.

USE: refrigerated or ambient temperature specimens.

** Styrofoam shippers measuring approx. 9" x 8" x 11" with a corrugated cardboard outer box.

USE: larger volumes of refrigerated or ambient temperature specimens or frozen (dry ice) specimens (note dry ice shipments are regulated by Federal and International laws and require special packaging, labeling, and paperwork).

MPLN Notes

Date Filled